


High-Throughput CellSensor® Reporter Assays For Multi-Pathway Analysis

Michael Hancock, Jun Wang, Jennifer Wilkinson, Thomas Machleidt, and Kun Bi
Invitrogen Corporation • Discovery Sciences • 501 Charmany Drive, Madison, WI 53719 • USA

Abstract

Recent studies of cell signaling pathways revealed many genetic changes that contribute to various diseases. We have developed a portfolio of cell-based reporter assays that can be used to screen for perturbagens of over 20 different endogenous signaling pathways as well as specific protein targets (e.g. kinases) that are involved in these pathways. To build these pathway reporter assays, we stably engineered the beta-lactamase reporter under the control of pathway-specific response elements into various cell backgrounds that natively express either wild-type or mutated signaling pathways. Activation of a pathway either through ligand binding or a clinically relevant constitutively active mutation in the pathway leads to the expression of beta-lactamase, which can be easily detected and quantified by a live cell, ratiometric beta-lactamase substrate. These cell lines were functionally validated with small molecule inhibitors and/or RNAi against known players of each pathway and were shown to be robust enough for high through-put screening. We have profiled 45 commercially available small molecule inhibitors against 11 signaling pathways. The results highlight the utility of multi-pathway cell-based profiling to advance drug discovery by providing robust data for both on and off pathway activities.

Figure 1 – GeneBLAZer® Beta-Lactamase Reporter Assay Technology



GeneBLAZer® Technology uses a FRET-based substrate to provide reliable and sensitive detection of beta-lactamase in cells. The advantages of this technology include:

1. Ratiometric readout reduces experimental noise
2. Fluorescent readout allows for live cell sorting by flow cytometry & imaging by microscopy
3. Sensitive detection due to enzyme turnover of substrate

Figure 2 – Invitrogen CellSensor® Cell Lines Offer Broad Pathway Coverage




Figure 3 – CellSensor® Cell Lines for Toxicity Pathway Analysis




Figure 4 – CellSensor® Cell Lines for Step-Wise Interrogation of Wnt/ β -catenin Pathway




Figure 5 – CellSensor® Cell Lines for Various Cytokine-mediated JAK/STAT Pathways




Figure 6 – CellSensor® Cell Lines for Compound Selectivity Profiling

IC50 (nM)	Pathway	Hypoxia	IFN α /JAK/STAT	IL-6/JAK/STAT	IL-1/IRAK	TNF α /NFKB	TGF β /SMAD	Wnt/PCSK	EGFR	IFN α /JAK/STAT	TGF β	PKC/Ca	cAMP/PKA
1000	IL-2	+	+	+	+	+	+	+	+	+	+	+	+
1000	IL-3	+	+	+	+	+	+	+	+	+	+	+	+
1000	IL-4	+	+	+	+	+	+	+	+	+	+	+	+
1000	IL-6	+	+	+	+	+	+	+	+	+	+	+	+
1000	GM-CSF	+	+	+	+	+	+	+	+	+	+	+	+
1000	EPO	+	+	+	+	+	+	+	+	+	+	+	+
1000	IFN α	+	+	+	+	+	+	+	+	+	+	+	+
1000	IFN β	+	+	+	+	+	+	+	+	+	+	+	+
1000	IFN γ	+	+	+	+	+	+	+	+	+	+	+	+
1000	Wnt3a	+	+	+	+	+	+	+	+	+	+	+	+
1000	EGFR	+	+	+	+	+	+	+	+	+	+	+	+
1000	PKC	+	+	+	+	+	+	+	+	+	+	+	+
1000	cAMP	+	+	+	+	+	+	+	+	+	+	+	+

Summary

- We have developed a panel of validation CellSensor® cell lines that can be used for endogenous signaling pathway analysis.
- Cell lines were further validated with specific Stealth™ RNAs or small molecular inhibitors against an essential player of the pathway.
- Use of CellSensor® cell lines to examine mutant and wild-type signaling pathways or to step-wise interrogate a pathway can facilitate a deeper understanding of disease biology.
- Many of our CellSensor® cell lines can be used for studying pathways initiated by various cytokines.
- CellSensor® cell lines can be used to generate compound potency and selectivity profiles.