


Human Adipose-Derived Stem Cells (ADSC)



Human Adipose-Derived Stem Cells (ADSC) were transfected using the Neon™ transfection system and 0.5 µg of a plasmid encoding the EGFP. 48 hours post transfection, the cells were analyzed by light (A) and fluorescence microscopy (B).

Electroporation parameters

Pulse voltage (v)	Pulse width (ms)	Pulse number	Cell density (cells/ml)	Transfection efficiency	Viability	Tip type
1,400	10	3	1 x 10 ⁷	88%	96%	10 µl
1,200	20	2	1 x 10 ⁷	79%	99%	10 µl

Cell information

Cell type/description	Human Adipose-Derived Stem Cells (ADSC)
Characteristics/species	Adherent/human
Tissue origin	Lipoaspirate
Media	Complete MesenPRO RS™ Medium with MesenPRO RS™ Growth Supplement and L-glutamine
Morphology	Spherical, colony
Double time	—
Subculturing	Medium renewal: every 2 to 3 days
Culture condition	Temperature: 37°C; atmosphere: air, 95%; carbon dioxide (CO ₂), 5%