

Comparison of Trypsin Solutions								
	TrypLE™ Select, Stable Trypsin-Replacement Enzyme (1X)	TrypLE™ Express, Stable Trypsin Replacement Enzyme (1X)	TrypLE™ Express, Stable Trypsin Replacement Enzyme (1X)	Trypsin, 0.25% (1X)	Trypsin, 2.5% (10X)	Trypsin-EDTA (0.05% Trypsin, EDTA•4Na) (1X)	Trypsin-EDTA (0.5% Trypsin, EDTA•4Na) (10X)	Trypsin-EDTA (0.25% Trypsin, EDTA•4Na) (1X)
	12563	12605	12604	15050	15090	25300	15400	25200
0.5 g/L of trypsin						●		
2.5 g/L of trypsin				●				●
5.0 g/L of trypsin							●	
25 g/L of trypsin					●			
1 M EDTA	●	●	●					
0.2 g/L of EDTA						●		
2.0 g/L of EDTA							●	
0.38 g/L of EDTA								●
w/ phenol red		●				●		●
w/o phenol red	●		●	●	●		●	
10X Concentration					●		●	
Animal Origin Free	●							