

Dedicated to the difference

AllSet⁺™ Gold SSP Typing Kits

- Dried, predispensed primer mixes for easy storage and use
- Predispensed amplification buffer saves time
- Negative contamination well for each test ensures quality
- UniMatch® Plus software and SCORE™ file aid in interpretation

Looking for an SSP typing alternative? *AllSet*⁺™ Gold SSP Kits are your answer.


These dry-format kits are simple to use, providing a high level of resolution and fewer ambiguities.

The *AllSet*⁺™ Gold SSP Kits are provided in precut polypropylene trays and come in a variety of kit formats:

- Direct to high-resolution kits for allele-level typing in a single tray (improve turnaround time for results, conserve on DNA samples, and reduce the number of kits to QC)
- Class I low-resolution kits, including combination trays
- Class I high-resolution kits
- Class II high-resolution kits

A free copy of UniMatch® Plus software is available with all orders.

www.invitrogen.com/transplantdiagnostics



Ordering information

Product	Quantity	Cat. no.
AllSet™ Gold SSP Typing Kits		
HLA-A High Resolution SSP	10 tests	54010D
HLA-B High Resolution SSP	10 tests	54020D
HLA-Cw High Resolution SSP	10 tests	54030D
HLA-DRB1 High Resolution SSP	10 tests	54040D
HLA-DQB1 High Resolution SSP	10 tests	54050D
HLA-A*11 High Resolution SSP	10 tests	541100D
HLA-A*02 High Resolution SSP	10 tests	54110D
HLA-A10 High Resolution SSP	10 tests	541110D
HLA-B17 High Resolution SSP	10 tests	541120D
HLA-A*03 High Resolution SSP	10 tests	541130D
HLA-A*01/36 SSP	10 tests	541140D
HLA-B*38/39 SSP	10 tests	541150D
HLA-B27 High Resolution SSP	10 tests	541160D
HLA-B22 SSP	10 tests	541170D
HLA-B*08 High Resolution SSP	10 tests	541180D
HLA-B*18 High Resolution SSP	10 tests	541190D
HLA-B*14 High Resolution SSP	10 tests	541200D
HLA-B*15 High Resolution SSP	10 tests	54120D
HLA-A*68 High Resolution SSP	10 tests	541210D
HLA-B*13 High Resolution SSP	10 tests	541220D
HLA-B21 High Resolution SSP	10 tests	541230D
HLA-A*19 High Resolution SSP	10 tests	54130D
HLA-B*35 High Resolution SSP	10 tests	54140D
HLA-B*44/45 High Resolution SSP	10 tests	54150D
HLA-A*23/24 SSP	10 tests	54160D
HLA-B*07 High Resolution SSP	10 tests	54170D
HLA-B*51/52 SSP	10 tests	54180D
HLA-B*40 High Resolution SSP	10 tests	54190D
KIR Genotyping	10 tests	54410D
HLA-A Low Resolution SSP	20 tests	54310D
HLA-B Low Resolution SSP	20 tests	54320D
HLA-ABC Low Resolution SSP	10 tests	54340D
HLA-Cw Low Resolution SSP	20 tests	54330D
HLA-ABDR Low Resolution SSP	10 tests	54350D
HLA-ABDRDQ Low Resolution SSP	10 tests	54360D
HLA-DR Low Resolution SSP	20 tests	54370D
HLA-DRDQ Low Resolution SSP	20 tests	54380D

Product	Quantity	Cat. no.
Molecular accessories		
Electro-Fast® gel accessories		
Electro-Fast® Gel System	1	920001
Fluorescent Numbering Panel	1	920011
Single 9-Tooth Comb	1	920021
End Dams	1 set	920031
Centipede™ gel accessories		
OWL Centipede™ Gel System	1	800001
Single 1.5 mm Comb	1	810001
Fluorescent Numbering Panel	1	810011
Casting Chamber	1	820011
Replacement Gaskets	1 pair	820021
Gasketed UV-Transparent Gel Tray	1	830001
Isolation kits		
DNA Isolation Kit	50 isolations	761001
Reagents		
PCR Markers, 50–2,000 bp	250 µl	74601250
PCR Markers, 100 bp ladder	250 µl	74602250
PCR Marker 50–500 bp	250 µl	74605250
Agarose	500 g	75000500
PCR Amplification Buffer GL - Plus	580 µl	741001
Heat Block	1	900001



www.invitrogen.com