

Scope of Services

The Invitrogen business of Life Technologies Corporation (“Invitrogen”) will profile a specified number of kinase(s) of interest in an *in vitro* Kinase Substrate Identification (KSI) assay using ProtoArray® microarrays that contain more than 9,000 human proteins as outlined below on behalf of the “Client”.

This Statement of Work specifies Invitrogen’s ProtoArray® KSI assay for this project.

Service includes:

1. Characterization of kinase substrates on ProtoArray® Human Protein Microarrays v5.0 (containing more than 9,000 unique proteins spotted in duplicate) in an *in vitro* kinase assay.
2. Invitrogen kinase, if chosen by the Client.
3. Controls and samples as listed below—this service includes six (6) ProtoArray® microarrays in total.
4. Probing of ProtoArrays® v5.0 at two (2) kinase concentrations (5 nM and 50 nM, unless otherwise specified by Client), in duplicate, and one (1) concentration of ³³P-ATP (33 nM).
5. Probing of one (1) ProtoArray® v5.0 with ³³P-ATP only as a negative control.
6. Probing of one (1) ProtoArray® v5.0 with MAPK14/p38α (100 nM) and ³³P-ATP.

Notes:

- a. *When multiple kinases are tested simultaneously, only one positive and negative control will be performed.*
- b. *Invitrogen does not guarantee that any interactions of the kinase of interest will be observed on the array or that observed interactions are specific to the kinase of interest.*

Service does not include:

1. Shipment of a ProtoArray® microarray.
2. Kinase with the exception of kinases available through the Invitrogen catalog.
3. Inhibitor compounds.
4. Assay using inhibitor compounds.
5. Assay performed at the same kinase activity for multiple kinases.

Timeline: 4 weeks upon initiation of service project

Client provides:

1. Kinase (at least two aliquots of 20 μl at 1 μM) or selection of kinase(s) from Invitrogen’s Catalog.
2. Signed Quotation and related Purchase Order.

Invitrogen delivers:

1. Comprehensive project report that includes all data of candidate kinase substrates, and results of the positive and negative controls for the ProtoArray® Kinase Substrate Identification Assay.
2. Image files and raw data.



ProtoArray[®] Kinase Substrate Profiling Assay

Ordering information:

Service Description	Catalog Number
ProtoArray [®] Kinase Substrate Profiling Service	Custom Quote

Timeline:

Upon receipt of all of the following: (i) this Quotation signed by Client; (ii) the Purchase Order related to this Quotation; and (iii) the Client Materials, Invitrogen will place the project in the current queue for execution. The Client will be informed of the predicted start date. Once the project is initiated, Invitrogen will use commercially reasonable efforts to supply the Deliverables to Client within the specified timelines. The timeline may be extended for orders of >3 kinases.

Project management:

Our project management team will be overseeing the work on this project and will keep you apprised of the status.

Terms and conditions:

Invitrogen's Terms and Conditions for Provision of Services govern the service work provided hereunder and are available upon request. Furthermore, this service may be covered by one or more Limited Use Label Licenses (See Technical Support section of Invitrogen.com for specific Limited Use Label Licenses). By use of this service, you accept the terms and conditions of all applicable Limited Use Label Licenses. This service is performed for research use only.

Contact information:
Invitrogen Custom Services Invitrogen, part of Life Technologies Corporation Phone: 800-955-6288 (3) x 45682 Fax: 716-774-3157 Email: custom.services@invitrogen.com